

# Technical specifications

Cycle description	Thermal disinfection temperature	Washing time	Thermal disinfection time	A <sub>0</sub>	Total cycle time
Handpiece cycle	90	13'	4'	3,000	34'



### Standard equipment

- 4 turbine adapters
- 3 micromotor handpiece adapters
- 1 scaler tip adapter (3 tips per adapter)
- 4 plugs

02/2021 SHMDGB211S00

Due to our policy of constant technological upgrading, technical specifications may be subject to change without prior notice. According to the standards in force, in extra-EU areas the availability and specifications of some products and/or characteristics may vary. Please contact your local distributor for further information.



## TETHYS HMD

### Revolutionary Evolution



Making Your Life Better.

**BU Medical Equipment**  
**Sede legale ed amministrativa**  
**Headquarters**  
 Cefla s.c.  
 Via Selice Provinciale, 23/a  
 40026 Imola - Bo (Italy)  
 tel. +39 0542 653111  
 fax +39 0542 653344

**Stabilimento**  
**Plant**  
 Via Bicocca, 14/c  
 40026 Imola - Bo (Italy)  
 tel. +39 0542 653441  
 fax +39 0542 653601



Making Your Life Better.





# Thanks to the new **HMD**, Tethys H10 is now also a handpiece thermal disinfecter

This latest accessory makes it possible to wash, thermally disinfect and dry up to 6 instruments (turbines or micromotor handpieces) simultaneously. The HMD is compatible with all Tethys H10, including those that have already been installed. All the user has to do is update the firmware to make the thermal disinfecter more complete than ever.



### Handpiece cycle

With the new handpiece cycle and the dedicated support, a cycle can be completed in just 34 minutes. This cycle consists of:

- Washing (no ultrasound)
- Thermal disinfection
- Drying

Thermal disinfection effectiveness on external surfaces has been verified as per EN15883-1 and EN15883-2.



### Maximum flexibility

The support can easily be assembled and removed; this allows use of a traditional basket and ensures the Tethys H10 can still perform its original function as an ultrasound thermal disinfecter.



### Adapters for every requirement

Users can treat the following instruments:

- Turbines
- Micromotor handpieces
- Scaler tips

To ensure proper instrument treatment all you need to do is fit the specific adapter for the relative instrument family.



**Micromotor handpiece adapter**  
Allows micromotor handpieces to be treated



**Turbine adapter**  
Allows turbines to be treated



**Scaler tip adapter**  
Allows simultaneous treatment of up to 3 scaler tips



**Plug**  
Necessary to seal unused connectors

