

# KaVo MASTERmatic LUX M45 L

Puts your ergonomics and patients' safety first.

Achieve best results, even in hard-to reach areas thanks to the unique head/knee angle combination of the new KaVo MASTERmatic LUX M45 L.





# The new addition to KaVo's family of reds

## Flexible treatment options

Compatibility with drill lengths of up to 25 mm.

## Unique head/knee angle combination compared to others

For great access and visibility in the molar area.



## Optimised 3-nozzle spray

Wets the bur along the entire length and ensures particularly homogenous cooling effect.

## Shedding some light on the subject

The tried and tested KaVo glass fiber optic technology illuminates every treatment area down to the smallest detail.

## Too cool to overheat

The KaVo CoolHead technology ensures a cool instrument head at all times.



	MASTERmatic LUX M45 L	MASTERmatic LUX M25 L	MASTERmatic LUX M05 L mini	EXPERTmatic LUX E25 L	EXPERTmatic E25 C
Material no.	1.014.0280	1.009.3630	1.009.3640	1.007.5550	1.007.5551
Angulation	120°/19°	100°/19°	100°/19°	90°/19°	90°/19°
Transmission ratio	1:4.2	1:5	1:5	1:5	1:5
Speed	168,000 rpm	200,000 rpm	200,000 rpm	200,000 rpm	200,000 rpm
Glass rod light conductor	with up to 25,000 LUX	with up to 25,000 LUX	with up to 25,000 LUX	with up to 25,000 LUX	-